Dental-Zylinderampullenspritze

REF 2666-18 / 2666-18US / 2666-22US / 2668-18 / 2668-18US

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
Maintenance Instructions

Dear Customer!

With the purchase of SCHWERT Dental Cartridge Syringes you have chosen high quality products, the correct handling and use of which are described in the following. In order to keep risks and unnecessary burdening of the patient, the user and third parties as low as possible, we request that you carefully look through these instructions for use and keep them for further reference.

Procedure

AUTOMATIC REPROCESSING

Products

Dental cartridge syringes; syringes for intraligamental anesthetics; Application syringe for elastic impression materials.

Instruction

Due to the constructive design of medical devices, a clearly defined limit of expected reprocessing cycles can not be determined. The service life of a medical device depends on the appropriate utilization of the product and its careful handling and maintenance.

Prior to reshipment of products for repair at the manufacturer's site, medical devices must be cleaned as well as sterilized according to the prescribed reprocessing requirements.

Preparations at the place of installation

Remove obvious dirt particles from medical devices immediately after utilization. It is not allowed to use fused agents or hot water (> 40°C) for cleaning as these agents lock residues and impede an effective purification process.

Transport

Secure storage as well as secure transport in closed boxes is recommended in order to avoid any damage to the medical device as well as environmental contamination.

Preparations for decontamination

For further reprocessing, medical devices must be disassembled or opened up as far as possible. For details see separate instructions.

Manual pre-cleaning

The medical device must be disassembled (see above) and cleaned by means of a soft multi-purpose brush and cold tap water until all obvious residues and defilements are removed. Lumen, drill holes and thread turns must be flushed for at least 10 seconds at a pressure of 3.8 bars by means of a water pistol.

Cleaning

Disassembled or opened-up medical devices must be put in a screen basin on the plug-in mode. Then start the cleaning process.

1. 4 minutes pre-washing by means of cold water
2. Draining
3. 5 minutes pre-washing at 55°C with 0.5% Neodisher Mediclean, Dr. Weigert (Hamburg)
4. Draining
5. 3 minutes neutralization with warm tap water (> 40°C)
6. Draining
7. 2 minutes intermediate flushing by means of warm tap water (> 40°C)
8. Draining

Desinfection

The automatic thermal disinfection must be carried out according to the national requirements regarding the A0-coefficient (see ISO 15883).

Drying

Drying of the exterior parts of medical devices by means of the drying cycle of the cleaning / disinfection equipment. Additional manual drying is possible by means of a lint-free cloth. Hollow space of medical devices must be dried by means of sterile compressed air.

Functional test, maintenance

Afterwards, a visual appraisal regarding purity is necessary. The assembling of medical devices, maintenance and functional tests must be carried out according to the instruction manual.

If necessary, the reprocessing procedure must be re-run until the medical device is obviously clean.

SCHWERT
SWORD · ESPADA
EPEE · SPADA

Dental Cartridge Syringes

REF 2666-18 / 2666-18US / 2666-22US / 2668-18 / 2668-18US

Packaging

Packaging of medical devices for sterilization must comply with standards according to ISO 11607 and EN 868.

Sterilization (autoclaving)

Sterilization of products by means of the fractionated pre-vacuum method (according to ISO 13060 / ISO 17665) in compliance with the individual national requirements.

1. Fractionated pre-vacuum method (3-fold)
2. Sterilization temperature of 134°C
3. Minimum exposure time: 3 minutes (full cycle)
4. Drying time: at least 10 minutes

Storage

Storage of sterilized medical devices must be carried out in a dry, clean and dust-free environment at a moderate temperature of 5°C - 40°C.

Information on validation of processing

The following test instructions, materials and machines have been used for validation:

Cleaning supplies:
Neodisher Mediclean (alkaline); Dr. Weigert (Hamburg)

Cleaning / disinfection equipment:
Miele G 7735 CD with plug-in module, vario-TD-program (without disinfection)
For details see report Cleaning: 17607011411 - 1
Sterilization: 17607010811 - 1

Additional instructions

According to the Ordinance on Installation, Operation and use of Medical Devices (MPBetreibV in German), the operator is responsible for validating any reprocessing procedures, even if the above described chemicals and machines are not available.

The operator must ensure that any reprocessing procedures – including resources, material and staff – are capable of achieving the required results. National legislation as well as the state of technology requires compliance with validated procedures.

Guarantee

A. Schweickhardt GmbH & Co. KG supplies only tested and defect-free products to its customers. All of our products are designed and manufactured to comply with the most stringent quality requirements. We assume no liability for products that have been altered compared to the original, used for purposes for which they were not intended or used improperly. Products that are repaired or altered at any company other than SCHWERT or a company commissioned for this purpose by SCHWERT are no longer authorized to bear the CE symbol.


EVIDENCE OF DECONTAMINATION

On the basis of statutory regulations and above all to protect our employees, we need a signed "Evidence of decontamination" with every product return.

Please ensure that this "Evidence of decontamination" is enclosed filled out and signed with every product return (complaint / repair / other reason for return) and package the product such that there is no risk of injury for our personnel in goods receiving during unpacking.

The proof of decontamination form "080300008a_DEKONTAMINATIONSNAHWEIS_AS.pdf" is available for download in the media library on our website "www.schwert.com".

The above instructions have been validated as SUITABLE for preparation of a medical device for its reuse by the medical device manufacturer. The reprocessor bears responsibility for achieving the desired results of the actual reprocessing with the equipment, materials and personnel used in the reprocessing facility for that purpose. This normally requires validation and routine monitoring of the process. Also, any deviation from the instructions provided on the part of the reprocessor should be carefully assessed for effectiveness and possible negative effects.


2666-18

1.8 ccm
Metric thread
Without aspiration device


2666-18US / 2666-22US

1.8 ccm / 2,2 ccm
US-thread
Without aspiration device

2668-18

1.8 ccm
Metric thread
With aspiration device

2668-18US

1.8 ccm
Metric thread
With aspiration device